

Product Description SALSA® MLPA® probemix P361-A3/P362-A3 USH2A

To be used with the MLPA General Protocol.

P361 Version-A3. As compared to version A2, three reference probes have been replaced. For complete product history see page 9.

P362 Version-A3. As compared to version A2, two reference probes have been replaced. For complete product history see page 9.

Catalogue numbers:

- **P361-025R:** SALSA® MLPA® probemix P361 USH2A mix 1, 25 reactions.
- **P361-050R:** SALSA® MLPA® probemix P361 USH2A mix 1, 50 reactions.
- **P362-025R:** SALSA® MLPA® probemix P362 USH2A mix 2, 25 reactions.
- **P362-050R:** SALSA® MLPA® probemix P362 USH2A mix 2, 50 reactions.

To be used in combination with a SALSA® MLPA® reagent kit, available for various number of reactions. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman capillary sequencers, respectively (see www.mlpa.com).

Certificate of Analysis: Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mlpa.com.

Precautions and warnings: For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mlpa.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

General Information: The SALSA MLPA Probemix P361 and P362 USH2A is a **research use only (RUO)** assay for the detection of deletions or duplications in the *USH2A* gene.

Usher syndrome is a clinically and genetically heterogeneous autosomal recessive disorder characterized by sensorineural hearing deficiencies at birth and later development of progressive retinitis pigmentosa (RP). It is the most frequent cause of combined deafness and blindness in adults and affects 3 to 6% of children born with hearing impairment. In brief, patients with Usher syndrome type II have mild hearing impairment with normal vestibular responses. Type II is the most common form of three Usher syndromes. Mutations within the *USH2A* gene have been associated with Usher syndrome type IIa and retinitis pigmentosa.

More information is available at <https://www.ncbi.nlm.nih.gov/books/NBK1341>

This SALSA MLPA Probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.

Gene structure and Transcript variants:

Entrez Gene shows transcript variants of each gene: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene>

For NM_ mRNA reference sequences: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide>

Locus Reference Genomic (LRG) database: <http://www.lrg-sequence.org/>

Exon numbering: The exon numbering used in this P361-A3/P362-A3 USH2A product description is the exon numbering from the RefSeq transcript NM_206933.2 for *USH2A*. The exon numbering and NM sequence used have been retrieved in July 2019. As changes to the NCBI database can occur after release of this product description, exon numbering may not be up-to-date.

Probemix content: The SALSA MLPA probemixes P361-A3 USH2A mix 1 and P362-A3 USH2A mix 2 both contain 45 probes with amplification products between 130 nt and 454 nt, and between 130 nt and 477 nt respectively. The *USH2A* gene (72 exons) spans ~800 kb of genomic DNA, and is located on chromosome 1q41, ~216 Mb from the p-telomere. The P361-A3 and P362-A3 probemixes together contain one probe for

each exon of the *USH2A* gene. The P361-A3 probemix contains probes for the odd exon numbers of *USH2A* and the P362-A3 probemix contains probes for the even exon numbers of *USH2A*. In both probemixes, nine reference probes are included, that detect different autosomal chromosomal locations.

These probemixes contain nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity Fragments (Q-fragments), two DNA Denaturation Fragments (D-fragments), one benchmark fragment, one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mlpa.com.

Length (nt)	Name
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)
88-96	D-fragments (low signal of 88 nt and 96 nt fragment indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

MLPA technique: The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mlpa.com).

MLPA technique validation: Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤ 0.10 for all probes over the experiment.

Required specimens: Extracted DNA from peripheral blood, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples: A sufficient number (≥ 3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from unrelated individuals who are from families without a history of sensorineural hearing disease. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol.

Positive control DNA samples: MRC-Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Biobank (<https://catalog.coriell.org>) and DSMZ (<https://www.dsmz.de/home.html>) have a diverse collection of biological resources which may be used as a positive control DNA sample in your MLPA experiments. The quality of cell lines can change, therefore samples should be validated before use.

Data analysis: Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mlpa.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results: The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 and the dosage quotient (DQ) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the DQ of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status	Dosage quotient
Normal	$0.80 < DQ < 1.20$
Homozygous deletion	$DQ = 0$
Heterozygous deletion	$0.40 < DQ < 0.65$
Heterozygous duplication	$1.30 < DQ < 1.65$
Heterozygous triplication/Homozygous duplication	$1.75 < DQ < 2.15$
Ambiguous copy number	All other values

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Incomplete DNA denaturation (e.g. due to salt contamination) can lead to a decreased probe signal, in particular for probes located in or near a GC-rich region or in or near the *USH2A* gene. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- When running MLPA products, the capillary electrophoresis protocol may need optimization. False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: lower injection voltage / injection time settings, or a reduced amount of sample by diluting PCR products.

Limitations of the procedure:

- In most populations, the major cause of genetic defects in the *USH2A* gene are small (point) mutations, most of which will not be detected by using SALSA MLPA Probemix P361-A3 and P362-A3.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNPs, point mutations, small indels) in the target sequence detected by a probe can cause false positive results. Mutations/SNPs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

Confirmation of results: Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

NCL gene variant database:

<https://databases.lovd.nl/shared/genes/USH2A>. We strongly encourage users to deposit positive results in the Leiden Open Variation Database (LOVD). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report copy number changes detected by the reference probes, false positive results due to SNPs and unusual results (e.g., a duplication of *USH2A* exons 6 and 8 but not exon 7) to MRC-Holland: info@mlpa.com.

Table 1a. SALSA MLPA Probemix P361-A3 USH2A mix 1

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) ^a	
		Reference	USH2A
64-105	Control fragments – see table in probemix content section for more information		
130	Reference probe 00797-L13645	5q31	
136	USH2A probe 11385-L12111		Exon 5
142	USH2A probe 11386-L12112		Exon 57
148	USH2A probe 11387-L12113		Exon 31
154	USH2A probe 11388-L12114		Exon 25
160	Reference probe 02493-L20011	17q11	
166	USH2A probe 11389-L12115		Exon 9
173	USH2A probe 11390-L20627		Exon 71
178	USH2A probe 11391-L12117		Exon 13
184	USH2A probe 11392-L12118		Exon 67
190	Reference probe 09561-L19082	20p13	
196	USH2A probe 11393-L12119		Exon 23
202	USH2A probe 11394-L12120		Exon 29
208	USH2A probe 16841-L19625		Exon 43
214	USH2A probe 11396-L12122		Exon 53
220 *	Reference probe 21710-L30368	15q21	
226	USH2A probe 11397-L12123		Exon 1
232	USH2A probe 11398-L12124		Exon 51
238	USH2A probe 11399-L12125		Exon 27
244	USH2A probe 11400-L12126		Exon 3
250 *	Reference probe 17332-L24586	3p21	
256	USH2A probe 11401-L12127		Exon 37
261	USH2A probe 11402-L12128		Exon 19
268	USH2A probe 16840-L29088		Exon 41
274	USH2A probe 11403-L12129		Exon 61
283	USH2A probe 11404-L12130		Exon 11
292 *	Reference probe 17041-L20104	10p11	
301	USH2A probe 11405-L12131		Exon 39
311	USH2A probe 11406-L20628		Exon 7
319	USH2A probe 11407-L12133		Exon 21
328	USH2A probe 11408-L12134		Exon 55
342	USH2A probe 20574-L20015		Exon 47
348	Reference probe 10713-L29089	6p12	
355	USH2A probe 11410-L12136		Exon 69
366	USH2A probe 11411-L20629		Exon 15
374	USH2A probe 11412-L12138		Exon 65
382	USH2A probe 11409-L19614		Exon 17
391	USH2A probe 11413-L12139		Exon 33
400	USH2A probe 11414-L12140		Exon 63
409	USH2A probe 11415-L12141		Exon 45
418	USH2A probe 11416-L12142		Exon 59
427	Reference probe 14774-L16471	11q23	
436	USH2A probe 11417-L12143		Exon 35
445	USH2A probe 11418-L12144		Exon 49
454	Reference probe 10061-L10485	8q22	

a) See above section on exon numbering for more information.

* New in version A3.

Table 1b. SALSA MLPA Probemix P362-A3 USH2A mix 2

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) ^a	
		Reference	USH2A
64-105	Control fragments – see table in probemix content section for more information		
130	Reference probe 00797-L13645	5q31	
136	USH2A probe 11421-L12147		Exon 50
142	USH2A probe 11422-L12148		Exon 40
148	USH2A probe 16839-L19623		Exon 32
154	USH2A probe 11424-L12150		Exon 60
160	Reference probe 12675-L13750	8q22	
166	USH2A probe 11425-L12151		Exon 18
172	USH2A probe 11426-L12152		Exon 30
178	USH2A probe 11427-L12153		Exon 62
184	USH2A probe 11428-L12154		Exon 14
190	USH2A probe 20575-L12157		Exon 64
196	USH2A probe 11429-L12155		Exon 66
202	USH2A probe 11430-L12156		Exon 42
214	USH2A probe 11432-L12158		Exon 48
220	Reference probe 16552-L19043	11q13	
226	USH2A probe 11433-L12160		Exon 70
232	USH2A probe 11434-L20182		Exon 54
238	USH2A probe 11435-L12163		Exon 68
244	USH2A probe 11436-L12165		Exon 4
250 *	Reference probe 17332-L24586	3p21	
256	USH2A probe 11437-L12167		Exon 24
265	USH2A probe 11438-L12168		Exon 6
274	USH2A probe 16838-L19622		Exon 22
283 *	Reference probe 13577-L23176	19p13	
292	USH2A probe 11440-L19615		Exon 36
301	USH2A probe 11441-L12171		Exon 46
310	USH2A probe 16836-L19620		Exon 12
319	USH2A probe 20576-L19618		Exon 72
330	USH2A probe 20577-L12173		Exon 38
337	USH2A probe 11444-L19619		Exon 2
346	Reference probe 10713-L11295	6p12	
358	USH2A probe 20578-L19621		Exon 20
364	USH2A probe 16843-L29091		Exon 52
373	USH2A probe 11448-L12178		Exon 10
382	Reference probe 13329-L14755	18q21	
391	USH2A probe 11449-L12179		Exon 44
400	USH2A probe 16835-L19617		Exon 58
409	USH2A probe 11451-L12181		Exon 28
419	USH2A probe 11452-L12182		Exon 8
427	Reference probe 13057-L14240	15q14	
436	USH2A probe 11453-L12183		Exon 56
447	USH2A probe 20579-L12184		Exon 16
454	USH2A probe 11455-L12185		Exon 34
463	USH2A probe 11456-L12186		Exon 26
477	Reference probe 20580-L14666	9q21	

a) See above section on exon numbering for more information.

* New in version A3.

Table 2. USH2A probes arranged according to chromosomal location

Length (nt) P361 / P362	SALSA MLPA probe	USH2A exon ^a	Ligation site NM_206933.2	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		<i>Start Codon</i>	<i>388-390 (Exon 2)</i>		
226	11397-L12123	Exon 1	110-111	CGTTGCCAGGCA-ACGAGGGACAGC	0.8 kb
337	11444-L19619	Exon 2	213-214	GTCTTCGCCATT-GCTAAGAACGTC	4.0 kb
244	11400-L12126	Exon 3	1000-1001	TGACACTGGGGA-GAATTCTTGTGA	53.6 kb
244	11436-L12165	Exon 4	1146-1147	GCATCTGGTACT-GTGCAAATAGGA	37.3 kb
136	11385-L12111	Exon 5	1186-1187	TAGAGCAGTTTG-TCGGAAGAATGC	2.2 kb
265	11438-L12168	Exon 6	1355-1356	CATTCTAATGA-TGCAGGAGACAC	1.3 kb
311	11406-L20628	Exon 7	1666-1667	GAATGAAAAACA-ATGGAGATTTGG	0.5 kb
419	11452-L12182	Exon 8	1748-1749	CAATGTCACATT-TAGCATCTGAC	1.7 kb
166	11389-L12115	Exon 9	1963-1964	GTCATGCCGATA-ATCTCGACACAA	29.7 kb
373	11448-L12178	Exon 10	2179-2180	CTTTTGAGCACT-TCAGAGGGGGAG	2.8 kb
283	11404-L12130	Exon 11	2267-2268	CTTTTTCCGACA-AGTTGGTGCAGA	38.3 kb
310	16836-L19620	Exon 12	2409-2410	AGGCAGTGCAAT-CAGTGCCAGAAT	3.9 kb
178	11391-L12117	Exon 13	2649-2650	CTCCATGGCTCA-GTGAACAAATTC	15.1 kb
184	11428-L12154	Exon 14	3308-3309	TGTTTGCCAAGA-TGCTTCCATTGC	14.6 kb
366	11411-L20629	Exon 15	3491-3492	GTGTGATGCTTG-TGTTCCAGTGC	10.1 kb
447	20579-L12184	Exon 16	3594-3595	AGTTCTTCTGCT-ATCAATCTCTCC	7.4 kb
382	11409-L19614	Exon 17	3887-3888	CTCTGTGACACT-TACCTGGACAAC	1.5 kb
166	11425-L12151	Exon 18	4372-4373	CCATCACTGGCT-TGGAGCCATACA	1.7 kb
261	11402-L12128	Exon 19	4494-4495	ATCCCTCCTTCA-GTCTTTCCCTC	6.4 kb
358	20578-L19621	Exon 20	4659-4660	ACTGCTAAATCC-CAAGAACTATCT	15.0 kb
319	11407-L12133	Exon 21	4913-4914	AGAGTCATCTCT-ACCAGCTCTGAT	78.2 kb
274	16838-L19622	Exon 22	5065-5066	TGATTGTCTTTG-CAGCATCACCTG	8.1 kb
196	11393-L12119	Exon 23	5231-5232	TGCTATTAGGCA-TCAGGCTTTTGG	2.3 kb
256	11437-L12167	Exon 24	5301-5302	AATGGTAGTACT-GTTATTGGAGAT	2.0 kb
154	11388-L12114	Exon 25	5488-5489	AAGAACAATCA-ACGTGTATAACA	1.2 kb
463	11456-L12186	Exon 26	5588-5589	TATGTTTCATGG-TGGAATGAACTT	5.3 kb
238	11399-L12125	Exon 27	5852-5853	ACTGATGAAGCA-TGCATCGGAGTC	5.0 kb
409	11451-L12181	Exon 28	6053-6054	AGTCAATCTGGA-TGGATGCCTATC	0.3 kb
202	11394-L12120	Exon 29	6191-6192	GATAGCCTCGCA-TGAAGGAGGTTTC	2.7 kb
172	11426-L12152	Exon 30	6315-6316	GTGACCTGGGAT-GAACCTGTTGTC	21.6 kb
148	11387-L12113	Exon 31	6511-6512	CTGGCTGTACTG-AGAGCTCACATG	2.1 kb
148	16839-L19623	Exon 32	6690-6691	TATTCAGGCAGT-GAGGAGAACTAC	46.0 kb
391	11413-L12139	Exon 33	6816-6817	CAGCTGCCACCA-GAACACGTGGAT	1.6 kb
454	11455-L12185	Exon 34	7044-7045	ATCAAGCTGGGA-GTAAGTTTATTT	5.8 kb
436	11417-L12143	Exon 35	7174-7175	TCTCCTGGACTG-AGCCTGAATATC	22.3 kb
292	11440-L19615	Exon 36	7281-7282	GGATTTGCTCCT-TGGAGTTTACAT	5.3 kb
256	11401-L12127	Exon 37	7464-7465	AATGGACTCTTA-ACACACTCAGTC	30.6 kb
330	20577-L12173	Exon 38	7571-7572	AGAGACAAACCT-TTGGGTGCTCAT	33.9 kb
301	11405-L12131	Exon 39	7761-7762	TGGTCTACACCA-GCTCGTAATAAC	0.7 kb
142	11422-L12148	Exon 40	7941-7942	GGATTTGGCAGT-GCACATAGTTCT	11.5 kb
268	16840-L29088	Exon 41	8429-8430	GAGTCATTCCAT-GAGGTTTATTGA	9.7 kb
202	11430-L12156	Exon 42	8826-8827	TTACCTACCTAT-GTTACCACTCAC	1.0 kb
208	16841-L19625	Exon 43	8979-8978 reverse	TTTGATGCAAGT-GGCTGCTGGATT	10.8 kb
391	11449-L12179	Exon 44	9153-9154	ACAACGTTAGCT-GGTCTTCCAGAG	21.2 kb
409	11415-L12141	Exon 45	9403-9404	GAGTCCACAGCA-TCAACAGTGCAG	1.4 kb
301	11441-L12171	Exon 46	9475-9476	CAGAGTTGTCA-TCATCAACAGTA	6.4 kb
342	20574-L20015	Exon 47	9690-9689 reverse	GTGGTAATTTGG-GTTCCATTGCTT	21.0 kb
214	11432-L12158	Exon 48	9889-9890	GTGATGAGCTCT-GCAAGGCAGTGA	3.3 kb
445	11418-L12144	Exon 49	10066-10067	GTGGTGCCGAA-TACAGGAGGCAC	14.8 kb
136	11421-L12147	Exon 50	10250-10251	GAGGCTTCATGA-TGGCCATGGCCA	8.8 kb
232	11398-L12124	Exon 51	10386-10387	AATATGTCAGAT-ACCATATGCTGC	3.5 kb
364	16843-L29091	Exon 52	10746-10747	ACCATTTCATACA-GGGAGTGTA AAC	3.8 kb

Length (nt) P361 / P362	SALSA MLPA probe	USH2A exon ^a	Ligation site NM_206933.2	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
214	11396-L12122	Exon 53	10825-10826	TTTCTGCCTGGA-ACAGCTATGGGC	0.8 kb
232	11434-L20182	Exon 54	11122-11123	GTGCCACCAGTA-GCAAGGTAAGAG	2.0 kb
328	11408-L12134	Exon 55	11164-11165	TTCCGGAGAGCA-TCCTGCCACCAA	13.3 kb
436	11453-L12183	Exon 56	11364-11365	TTCACCTTTACA-GCTTGTACATCT	7.0 kb
142	11386-L12112	Exon 57	11475-11476	ATCAATTCTACA-ACAGTGGGAATTA	1.1 kb
400	16835-L19617	Exon 58	11665-11666	GCAGCAGTGCTA-GTGATGATTACA	15.4 kb
418	11416-L12142	Exon 59	11829-11830	CTCAATGATGGA-AGTGTAAACACCT	1.8 kb
154	11424-L12150	Exon 60	12027-12028	GCACTGGGGTCA-GCTTGCATAGAG	13.4 kb
274	11403-L12129	Exon 61	12391-12392	ACCGTGTGGTCT-ACCAGGAGAGAC	47.9 kb
178	11427-L12153	Exon 62	12641-12642	GGCATTGCTACT-ACAGTGGTCAGA	4.8 kb
400	11414-L12140	Exon 63	12882-12883	ACTGTCCACTCT-GTGAAGTCCACC	4.3 kb
190	20575-L12157	Exon 64	14331-14332	GCTCCAGGAGGA-TTCCAGCCAAC	20.5 kb
374	11412-L12138	Exon 65	14636-14637	CTCTTCTACCCA-AGCAGTGGTCAA	2.0 kb
196	11429-L12155	Exon 66	14821-14822	GCACCTGCTTCA-ACTGTTGCAGCA	1.0 kb
184	11392-L12118	Exon 67	15004-15005	CTTGCCCTCTG-ACTCAGCCCTCC	7.0 kb
238	11435-L12163	Exon 68	15230-15231	GTCTGTGGTGTG-TGTGAAGTGGAG	1.5 kb
355	11410-L12136	Exon 69	15385-15386	TCTGCACGACTG-ACGAAGGAAGTG	4.6 kb
226	11433-L12160	Exon 70	15518-15519	GCTGTGGTTCAT-AGTGTAAATGGC	5.7 kb
173	11390-L20627	Exon 71	15780-15781	AGTCAAAACCAA-ACCAGCCTAAC	3.1 kb
319	20576-L19618	Exon 72	15965-15966	AGTGACTAAGGA-ACGCACCACATT	
		Stop Codon	15994-15996 (Exon 72)		

a) See above section on exon numbering for more information.

b) Only partial probe sequences are shown. Complete probe sequences are available at www.mlpa.com. Please notify us of any mistakes: info@mlpa.com.

References


- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

Selected publications using SALSA MLPA Probemix P361/P362 USH2A

- Dad S et al. (2015). Partial USH2A deletions contribute to Usher syndrome in Denmark. *Eur J Hum Genet.* Dec;23(12):1750. doi: 10.1038/ejhg.2015.131.
- Garcia-Garcia G et al. (2014). Novel deletions involving the USH2A gene in patients with Usher syndrome and retinitis pigmentosa. *Mol Vis.* Sep 25;20:1398-410. eCollection 2014.
- Steele-Stallard et al. (2013). Screening for duplications, deletions and a common intronic mutation detects 35% of second mutations in patients with USH2A monoallelic mutations on Sanger sequencing. *Orphanet J Rare Dis.* 2013; 8: 122.

P361 Product history	
<i>Version</i>	<i>Modification</i>
A3	Three reference probes have been replaced.
A2	One reference probe is replaced and the length of some probes have been adjusted.
A1	First release.
P362 Product history	
<i>Version</i>	<i>Modification</i>
A3	Two reference probes have been replaced.
A2	The length of several probes has been adjusted.
A1	First release.

Implemented changes in the product description	
<i>Version A3/A3-01– 15 July 2019 (02P)</i>	
<ul style="list-style-type: none"> - Product description restructured and adapted to a new template. - Product description adapted to a new product version. - Small changes in Table 1 and Table 2. 	
<i>Version 05 – 06 April 2016 (55)</i>	
<ul style="list-style-type: none"> - Product description adapted to new product versions (version numbers changed, lot numbers added, small changes in Table 1 and Table 2, new pictures included). - Various minor textual changes. - Related SALSA MLPA probemix added on page 1. 	

More information: www.mlpa.com; www.mlpa.eu	
	MRC-Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, the Netherlands
E-mail	info@mlpa.com (information & technical questions); order@mlpa.com (orders)
Phone	+31 888 657 200