

# Product Description

## SALSA® MLPA® Probemix P018-G2 SHOX

To be used with the MLPA General Protocol.

### Version G2

For complete product history see page 15.

### Catalogue numbers:

- **P018-025R:** SALSA MLPA Probemix P018 SHOX, 25 reactions.
- **P018-050R:** SALSA MLPA Probemix P018 SHOX, 50 reactions.
- **P018-100R:** SALSA MLPA Probemix P018 SHOX, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see [www.mrcholland.com](http://www.mrcholland.com)).

### Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at [www.mrcholland.com](http://www.mrcholland.com).

### Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: [www.mrcholland.com](http://www.mrcholland.com). It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

**Please note that copy number variation in healthy individuals has been reported for several genes outside of the SHOX region.**

### Intended purpose

The SALSA MLPA Probemix P018 SHOX is an in vitro diagnostic (IVD)<sup>1</sup> or research use only (RUO) semi-quantitative assay<sup>2</sup> for the detection of deletions or duplications in the human short stature homeobox (*SHOX*) gene and its regulatory regions on Xp22.33/Yp11.32 in genomic DNA isolated from human peripheral whole blood specimens or buccal swabs. P018 SHOX is intended to confirm a potential cause for disorders associated with short stature, including Leri-Weill dyschondrosteosis (LWD), Langer mesomelic dysplasia (LMD) or Idiopathic short stature (ISS).

Copy number variations (CNVs) detected with P018 SHOX should be confirmed with a different technique. In particular, CNVs detected by only a single probe always require confirmation by another method. In the majority of patients, defects in the *SHOX* gene region are deletions. However, point mutations can occur which will not be detected by MLPA. It is therefore recommended to use this assay in combination with sequence analysis.

Assay results are intended to be used in conjunction with other clinical and diagnostic findings, consistent with professional standards of practice, including confirmation by alternative methods, clinical genetic evaluation, and counselling, as appropriate. The results of this test should be interpreted by a clinical molecular geneticist or equivalent.

This device is not intended to be used for standalone diagnostic purposes, pre-implantation or prenatal testing, population screening, or for the detection of, or screening for, acquired or somatic genetic aberrations.

<sup>1</sup> Please note that this probemix is for in vitro diagnostic (IVD) use in the countries specified at the end of this product description. In all other countries, the product is for research use only (RUO).

<sup>2</sup> To be used in combination with a SALSA MLPA Reagent Kit and Coffalyser.Net analysis software.

## Clinical background

*SHOX* is located in the pseudoautosomal region 1 (PAR1) on the short arm of the X and Y chromosomes. Located upstream and downstream of *SHOX* are highly conserved non-coding elements (CNEs), some of which have been shown to be important *SHOX* enhancer sequences. Mutations in *SHOX* or its regulatory regions cause a range of disorders associated with short stature, including LWD, LMD, and ISS, as *SHOX* is a known transcription factor highly expressed in tissues responsible for bone development (Benito-Sanz et al. 2012b).

LWD is a dominant skeletal disorder characterised by short stature, mesomelic shortening of the limbs, and the characteristic Madelung deformity. LMD is a more severe form of LWD and is a result of mutations in both *SHOX* alleles (Bertorelli et al. 2007, Campos-Barros et al. 2007, Shears et al. 2002, Zinn et al. 2002). ISS classifies individuals with a height below the third centile in whom no identifiable disorder is present. Heterozygous mutations of *SHOX* and/or its regulatory elements are detected in approximately 60% of LWD patients and approximately 5-15% of ISS cases. Homozygous or compound heterozygous mutations of *SHOX* and/or its enhancers are detected in 75% of LMD patients (Benito-Sanz et al. 2006, Benito-Sanz et al. 2012a, Chen et al. 2009, Huber et al. 2006).

In individuals with a *SHOX* related disorder, 70-80% of *SHOX* mutations are whole gene deletions, 2-6% are partial deletions, and 20-25% are point mutations (Binder 2011, Caliebe et al. 2012). Duplications have also been reported in LWD and ISS patients (Benito-Sanz et al. 2011b). The P018 *SHOX* probemix can detect most deletions and duplications and therefore complements sequence analysis of *SHOX*.

More information is available on <http://www.ncbi.nlm.nih.gov/books/NBK1215/>.

## Gene structure

The *SHOX* gene spans 35 kilobases (kb) of the pseudoautosomal region 1 (PAR1) located on Xp22.33 / Yp11.32 and contains 7 exons. The *SHOX* LRG\_710 is available at [www.lrg-sequence.org](http://www.lrg-sequence.org) and is identical to GenBank NG\_009385.2.

## Transcript variants

For *SHOX*, two major transcript variants have been described: <http://www.ncbi.nlm.nih.gov/gene/6473>. *SHOX* transcript variant 1 (NM\_000451.4, 7934 nt, coding sequence 108-986) represents the longer transcript and encodes the longer active isoform (SHOXa). The ATG translation start codon is located in exon 1 (LRG exon 2b) and the stop codon is located in exon 5 (LRG exon 6). *SHOX* transcript variant 2 (NM\_006883.2, 1951 nt, coding sequence 692-1369) contains alternative 5' and 3' exons compared to transcript variant 1 and encodes a shorter isoform (SHOXb) with a different C-terminus than isoform SHOXa. Several regulatory sequences located outside of *SHOX* that affect *SHOX* transcription have been described (e.g. Benito-Sanz et al. 2012b, Durand et al. 2010, Fukami et al. 2006, Sandbacka et al. 2011).

## Exon numbering

The *SHOX* exon numbering used in this P018-G2 *SHOX* product description is the exon numbering from the LRG\_710 sequence. This exon numbering is different from the *SHOX* exon numbering in many articles where exon 7 is referred to as exon 6b. The exon numbering of the NM\_ sequence that was used for determining a probe's ligation site does not always correspond to the exon numbering obtained from the LRG sequences. As changes to the databases can occur after release of this product description, the NM\_ sequence and exon numbering may not be up-to-date.

## Probemix content

The SALSA MLPA Probemix P018-G2 *SHOX* contains 48 MLPA probes with amplification products between 124 and 504 nucleotides (nt). This includes 32 probes for the PAR1 region on chromosome Xp22 / Yp11, including at least one probe for each exon of *SHOX* transcript variant 1, and one probe for intron 6 only present in the *SHOXb* splice variant. Several probes are present for *SHOX* regulatory regions, located upstream and downstream of *SHOX*. Moreover, this probemix contains multiple flanking probes targeting the X chromosome outside the *SHOX* area: one probe detecting the area just before the *SHOX* upstream regulatory regions, five probes inside the PAR1 region but downstream of the *SHOX* area, and seven probes outside of PAR1. Flanking

probes can be used to characterise larger deletions/duplications and to distinguish *SHOX* deletions from a Turner syndrome karyotype. In addition, nine reference probes are included that detect autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online ([www.mrcholland.com](http://www.mrcholland.com)).

This probemix contains ten quality control fragments generating amplification products between 64 and 118 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and two chromosome Y-specific fragments (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at [www.mrcholland.com](http://www.mrcholland.com).

Length (nt)	Name
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)
88-96	D-fragments (low signal indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific) ( <i>AMOT</i> gene)
105	Y-fragment (Y chromosome specific) ( <i>UTY</i> gene)
118	Y-fragments (Y chromosome specific) ( <i>ZFY</i> gene)

#### MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol ([www.mrcholland.com](http://www.mrcholland.com)).

#### MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals of the same sex is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation  $\leq 0.10$  for all reference probes over the experiment.

#### Required specimens

Extracted DNA from peripheral blood or buccal swab, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

#### Reference samples

A sufficient number ( $\geq 3$ ) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different unrelated individuals who are from families without a history of short stature. **Importantly, all reference samples need to be of the same sex for correct data analysis.** To facilitate interpretation, it is also recommended to use reference and patient samples of the same sex. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol ([www.mrcholland.com](http://www.mrcholland.com)).

#### Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ (<https://www.dsmz.de/>) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. Sample ID numbers NA20212, NA20217, NA20218, and NA04626 from the Coriell Institute have been tested with this P018-G2 probemix at MRC Holland and can be used as a positive control samples. For details about genotype, affected probes and expected results, please see table below. The quality of cell lines can change; therefore samples should be validated before use.

Sample ID	Genotype	Affected probes	Expected final ratio
NA20212	Female sample: Heterozygous deletion (~0.9 Mb) of the <i>SHOX</i> gene and the upstream and downstream <i>SHOX</i> area.	09333-L10292 to 09335-L30792	0.5
NA20217	Male sample: Compound heterozygous deletion of <i>SHOX</i> and of <i>SHOX</i> downstream area.	01341-L20651 to 09338-L24247 13297-L24253 to 09335-L30792	0.5 0.5
NA20218	Female sample: Compound heterozygous deletions of upstream and downstream area of <i>SHOX</i> resulting in a homozygous deletion of the entire <i>SHOX</i> gene.	09333-L10292, 05642-L05096 to 14697-L24245  18889-L25087 to 09338-L24247	0.5  0
NA04626	Female sample: Heterozygous duplication of the <i>SHOX</i> gene and the upstream and downstream <i>SHOX</i> area.	09333-L10292 to 01156-L00659	1.5

### Performance characteristics

In individuals with a *SHOX* related disorder, 70-80% of all mutations are whole gene deletions and 2-6% are partial deletions, both of which can be detected by this MLPA probemix (Binder 2011, Caliebe et al. 2012). Duplications of *SHOX* have also been reported in LWD and ISS patients (Benito-Sanz et al. 2011b) and can be detected by this probemix. The analytical sensitivity and specificity for the detection of deletions/duplications in *SHOX* and its surrounding enhancer regions is very high and can be considered >99% (based on a 2006-2022 literature review).

Analytical performance can be compromised by: SNVs or other polymorphisms in the DNA target sequence, impurities in the DNA sample, incomplete DNA denaturation, the use of insufficient or too much sample DNA, the use of insufficient or unsuitable reference samples, problems with capillary electrophoresis or a poor data normalisation procedure and other technical errors. The MLPA General Protocol contains technical guidelines and information on data evaluation/normalisation.

### Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at [www.mrcholland.com](http://www.mrcholland.com). Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

**Note: Slope correction in samples with large deletions:** The slope correction algorithm in Coffalyser.Net has been optimised to give the best possible results in as many situations as possible. However, the slope correction algorithm may confuse a large deletion for sloping, which can lead to over- or under-correction. Since large deletions are often detected when using the P018 probemix, this issue may occur relatively frequently. Incorrectly applied slope correction can cause an FSLP warning in Coffalyser.Net or ambiguous results for multiple probes. If you suspect that slope correction was incorrectly applied we recommend contacting [info@mrcholland.com](mailto:info@mrcholland.com) for assistance.

### Interpretation of results

The expected results for *SHOX*-specific MLPA probes are allele copy numbers of 2 (normal), 0 (homozygous deletion), 1 (heterozygous deletion), 3 (heterozygous duplication), and 4 (heterozygous triplication/homozygous duplication). Please see Table 3 for examples of potential results.

The standard deviation of each individual reference probe over all the reference samples should be  $\leq 0.10$  and the final ratio (FR) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results when **reference samples of the same sex** have been used:

Copy number status		Final ratio (FR)
Pseudoautosomal sequences in males and females and X chromosome sequences in females	X chromosome sequences in males	
Normal	Normal	$0.80 < FR < 1.20$
Homozygous deletion	Deletion	FR = 0
Heterozygous deletion		$0.40 < FR < 0.65$
Heterozygous duplication		$1.30 < FR < 1.65$
Heterozygous triplication/homozygous duplication	Duplication	$1.75 < FR < 2.15$
Ambiguous copy number		All other values

Note: The term “dosage quotient”, used in older product description versions, has been replaced by “final ratio” to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

**Please note that the above mentioned final ratios can be affected in mosaic cases.**

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

P018 specific notes:

- Short stature and skeletal deformities due to *SHOX* defects are pseudoautosomal dominant disorders. Therefore, a heterozygous mutation of *SHOX* is expected to result in these phenotypes.
- A recurrent 47.5 kb deletion downstream of *SHOX* has been described by Benito-Sanz et al. (2012b). This 47.5 kb deletion is covered by three probes in this P018-G2 probemix (Table 2).



- Complete or partial duplications found within *SHOX* or its surrounding regulatory regions have been found in LWD and ISS patients (Benito-Sanz et al. 2011b).
- Breakpoints of partial *SHOX* deletions have been reported to frequently occur in intron 3 (Benito-Sanz et al. 2017).
- A partial *SHOX* deletion encompassing the last 20 nt of *SHOX* exon 3 and part of intron 3 has been described (Funari et al. 2019). This deletion will not be detected by P018-G2 *SHOX* as there is no probe present that targets this region.
- Deletion of the *SHOX* exon 6 and intron 6 probes, which are located downstream of the stop codon of transcript variant 1 (*SHOXa*), may not affect *SHOX* gene function.
- Please note that single exon deletions have a considerable chance of being a false positive result, either due to non-pathogenic copy number variants (Benito-Sanz et al. 2011a), mutations within the probe binding site (Barca-Tierno et al. 2011) or due to impurities in the DNA.
- Not all copy number changes detected by *SHOX* AREA probes will affect *SHOX* gene function. Analysis of family members may be required for correct interpretation of results.
- Large deletions and duplications interrupted by one or more probes with a normal copy number have been reported and might be the result of an inversion followed by a deletion/duplication or may indicate compound heterozygosity (Dupont et al. 2007). In these cases, parental evaluation can assist data interpretation.
- Flanking probes have been included in this probemix to help determine the extent of a deletion/duplication. Copy number changes detected by flanking probes only have been reported in healthy individuals and are unlikely to be related to short stature. However, several genes detected by the flanking probes have been associated with other disorders (Balasubramanian and Crowley Jr 2017, Mehta and Ebens 2021, Mullighan et al. 2009, Nguyen et al. 2022, Russell et al. 2009).

#### Limitations of the procedure

- The SALSA MLPA Probemix P018 *SHOX* will not detect point mutations in the *SHOX* gene, which are the second most common cause of genetic defects in *SHOX*.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

#### Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

#### *SHOX* mutation database

<https://databases.lovd.nl/shared/genes/SHOX>. We strongly encourage users to deposit positive results in the Leiden Open Variation Database (LOVD). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results (e.g., a duplication of *SHOX* exons 3 and 5 but not exon 4) to MRC Holland: [info@mrcholland.com](mailto:info@mrcholland.com).

**Table 1. SALSA MLPA Probemix P018-G2 SHOX**

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) <sup>a</sup>		
		Reference	SHOX region/PAR1	Outside PAR1
64-118	Control fragments – see table in probemix content section for more information			
124	Reference probe 15370-L13762	7q		
130	<b>SHOX AREA probe</b> 18885-L24430		CNE-3	
136	<b>SHOX AREA probe</b> 05642-L05096		6.2 kb after CNE2	
142 « ~	<b>IL3RA probe</b> 13597-L15055		Xp22.33 / Yp11.32	
148	<b>SHOX AREA probe</b> 05648-L06218		Downstream	
154	<b>SHOX AREA probe</b> 13821-L14642		CNE3	
160	Reference probe 04966-L04696	1p		
166	<b>SHOX probe</b> 01145-L00702		<b>Exon 1</b>	
172	<b>SHOX AREA probe</b> 18886-L24431		CNE4	
178	<b>SHOX AREA probe</b> 05649-L20176		Downstream	
185	<b>SHOX AREA probe</b> 06293-L20177		CNE9	
191	Reference probe 06057-L05512	4p		
199	<b>SHOX AREA probe</b> 13296-L20175		CNE5	
204	<b>SHOX probe</b> 01146-L06220		<b>Exon 2</b>	
211 « ~	<b>PPP2R3B probe</b> 09333-L10292		Xp22.33 / Yp11.32	
219	Reference probe 03247-L02684	13q		
226 ◊	<b>SHOX probe</b> 09336-L20178		<b>Exon 6</b>	
231	<b>SHOX probe</b> 09337-L00911		<b>Exon 6</b>	
238 ~ x	<b>ANOS1 probe</b> 06402-L09795			Xp22.31
245	<b>SHOX probe</b> 01147-L00802		<b>Exon 3</b>	
254 ~ x	<b>ARSF probe</b> 16846-L20647			Xp22.33
261	Reference probe 00587-L20649	18q		
266	<b>SHOX probe</b> 01341-L20651		Upstream	
274 ~ x	<b>FANCB probe</b> 03906-L03066			Xp22.2
283 ~ x	<b>NLGN4X probe</b> 05587-L04577			Xp22.31
290	<b>SHOX AREA probe</b> 06291-L06222		CNE9	
300 j	<b>SHOX probe</b> 01148-L15501		<b>Exon 4</b>	
310 ~	<b>ASMT probe</b> 01153-L00712		Xp22.33 / Yp11.31	
318 ^	<b>SHOX AREA probe</b> 05645-L05099		2 kb before CNE8	
328 ~ x	<b>PRKX probe</b> 16898-L19768			Xp22.33
337 ¥	<b>SHOX probe</b> 21538-L30066		<b>Exon 5</b>	
346	Reference probe 06560-L06118	1q		
355 ~ +	<b>VAMP7 probe</b> 01156-L00659			Xq28 / Yq12 (PAR2)
364	<b>SHOX AREA probe</b> 18889-L25087		CNE-5	
379	<b>SHOX AREA probe</b> 14697-L24245		Downstream	
389 ~	<b>CSF2RA probe</b> 10251-L24246		Xp22.33 / Yp11.32	
395	<b>SHOX probe</b> 09338-L24247		Intron 6	
403 ~	<b>CRLF2 probe</b> 13911-L19678		Xp22.33 / Yp11.32	
412	Reference probe 09793-L12593	15q		
420 ~ x	<b>AIFM1 probe</b> 00820-L25090			Xq25
427	<b>SHOX AREA probe</b> 18891-L25088		CNE-2	
439 ^	<b>SHOX AREA probe</b> 05646-L24249		5.4 kb after CNE8	
445 ¥	<b>SHOX AREA probe</b> 09335-L30792		Downstream	
456 ~	<b>ZBED1 probe</b> 16858-L25227		Xp22.33 / Yp11.31	
466	<b>SHOX AREA probe</b> 13297-L24253		6 kb before CNE9	
476	Reference probe 09888-L10301	16p		
481 ^	<b>SHOX AREA probe</b> 18893-L25091		CNE7	
504	Reference probe 09870-L19465	2p		

<sup>a</sup> See section Exon numbering on page 2 for more information.



**CNE** = Conserved Non-coding DNA Element. Locations of the upstream regulatory regions (CNE-2, CNE-3 and CNE-5) are based on Durand et al. (2010) and Benito-Sanz et al. (2012b). Locations of the downstream regulatory regions (CNE2 through CNE9) are based on Benito-Sanz et al. (2012b) and Fukami et al. (2006).

¥ Changed in version G2. Minor alteration, no change in sequence detected.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

- Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

x X-chromosome, outside PAR region. Gives 50% reduced signal in males as compared to females.

∫ A founder *SHOX* point mutation in the Spanish Gypsy population, c.508G>C (p.A170P), results in a strongly reduced probe signal of the 300 nt exon 4 probe (Barca-Tierno et al. 2011).

^ The 481, 318 and 439 nt probes are within the recurrent ~47.5 kb deletion region described by Benito-Sanz et al. (2012b).

+ The VAMP7 probe at 355 nt is located very close to the X and Y q-telomere in PAR2.

◇ From product description version 05 onwards the 226 nt probe is considered an exon 6 probe to align with the LRG exon numbering.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Please note: not all known SNVs are mentioned in the tables above. Single probe aberration(s) must be confirmed by another method.

**Table 2. Xp / Yp probes arranged according to chromosomal location**

Length (nt)	SALSA MLPA probe	Gene/Exon <sup>a</sup>	Ligation site	Partial sequence <sup>b</sup> (24 nt adjacent to ligation site)	Distance to next probe
<b>p-telomere</b>					
----- <b>Start of PAR1</b> -----					227.0 kb
211 « ~	09333-L10292	PPP2R3B gene	Xp22.33 / Yp11.32	CGTCCGAGTTCC-ACTCGCGCTACA	91.6 kb
<b>SHOX region</b>					
364	18889-L25087	SHOX AREA	CNE-5	GAAATGTTAACA-GCTCCCCGAGCT	61.7 kb
130	18885-L24430	SHOX AREA	CNE-3	ATGGCAGAGCAT-TTGTACCCCTGG	56.3 kb
427	18891-L25088	SHOX AREA	CNE-2	TACACCGTTATG-CGGATGCTCGTT	63.5 kb
266	01341-L20651	SHOX LOC159015	Upstream of SHOX	GCCTGGAACAGA-ACTTCCGCGGGG	4.7 kb
		<b>SHOX</b>	<b>NM_000451.4</b>		
		<i>start codon</i>	108-110 (Exon 2b)		
166	01145-L00702	Exon 1	6.3 kb before exon 2b (NM_006883.2: 99-100)	TTTCTACTGCAA-ACAGAAATGGGA	6.7 kb
204	01146-L06220	Exon 2	336-337	ACCACGTAGACA-ATGACAAGGAGA	3.6 kb
245	01147-L00802	Exon 3	448-449	CGGGCAGACCAA-GCTGAAACAGAG	6.2 kb
300 j	01148-L15501	Exon 4	614-615	CAGAACCGGAGA-GCCAAGTGCCGC	0.2 kb
337	21538-L30066	Exon 5	677-678	ACAGCCAACCAC-CTAGACGCCTGC	3.6 kb
231	09337-L00911	Exon 6	922-923	AAGCAACAGCAA-GAATTCAGCAT	6.4 kb
226 d	09336-L20178	Exon 6	7370-7371	TGGCTTCACGAG-TTCAGCCCATTG	6.4 kb
		<i>stop codon</i>	984-986 (Exon 6)		
395	09338-L24247	SHOX Intron 6	1.4 kb before exon 7	TCCCACATTCTT-GGAATCACAATG	56.9 kb
136	05642-L05096	SHOX AREA	6.2 kb after CNE2	GCAGCAGTGAAA-GTGAGCATTCCC	19.8 kb
154	13821-L14642	SHOX AREA	CNE3	GATGGCTGATAA-TTACTCCGTATG	19.4 kb
172	18886-L24431	SHOX AREA	CNE4	GCCTCCGATACA-GTTTACGGCTTC	37.4 kb
199	13296-L20175	SHOX AREA	CNE5	GGAAAACCACGT-TCCTATCGATCC	29.6 kb
481 ^	18893-L25091	SHOX AREA	CNE7	CAGACCAGGTCT-CCTGTTTCATGT	28.1 kb
318 ^	05645-L05099	SHOX AREA	2 kb before CNE8	TGTTCCACCGT-AAAACTCACTCC	8.4 kb
439 ^	05646-L24249	SHOX AREA	5.4 kb after CNE8	TGCATGTCTGCT-TTTTGAATGGCC	10.7 kb
466	13297-L24253	SHOX AREA	6 kb before CNE9	TACAGCAAATGA-TACGTATAAATT	6.3 kb
290	06291-L06222	SHOX AREA	CNE9	CTTGAAAGGGCA-GGAACTCTAATT	0.4 kb
185	06293-L20177	SHOX AREA	CNE9	TAATTGATGAGA-TGCAGAAGCCAG	15.4 kb
148	05648-L06218	Xp22-PAR1	Downstream	TGGTGCTGAAAT-GAGGAAGCCCTG	48.7 kb
178	05649-L20176	Xp22-PAR1	Downstream	TGAGGAGGTACC-TCAAAGCTAAAC	64.4 kb
445	09335-L30792	Xp22-PAR1	Downstream	GAAATTCAGTTT-TAATAACACAGA	66.0 kb
379	14697-L24245	Xp22-PAR1	Downstream	CTCTGGTGAGAT-GCCATCTAGAGA	338.0 kb
<b>End of PAR1</b>					
403 ~	13911-L19678	CRLF2 gene	Xp22.33 / Yp11.32	GAATGCCAGCAA-ATACTCCAGGAC	73.9 kb
389 ~	10251-L24246	CSF2RA gene	Xp22.33 / Yp11.32	GACAAGCCTTCT-GCTCTGTGAGTT	69.8 kb
142 « ~	13597-L15055	IL3RA gene	Xp22.33 / Yp11.32	TGCACAGATAAG-TTTGTCGTCTTT	280.7 kb
310 ~	01153-L00712	ASMT gene	Xp22.33 / Yp11.31	GACATCCAGAA-GTGGTGTGGACG	706.4 kb
456 ~	16858-L25227	ZBED1 gene	Xp22.33 / Yp11.31	TCGTCAAGAGCA-ACACGGAGCAGA	593.8 kb
254 ~ x	16846-L20647	ARSF gene	Xp22.33	CATCCATATAAT-TATGGGTTTGAC	536.9 kb
328 ~ x	16898-L19768	PRKX gene	Xp22.33	CGATTAGGAAAC-ATGAAGGTCAGT	2.6 Mb
283 ~ x #	05587-L04577	NLGN4X gene	Xp22.31	GACGGCTTGGGT-GATGCACGAAAT	2.3 Mb
238 ~ x #	06402-L09795	ANOS1 gene	Xp22.31	GTTTCCTGAAGC-GTGTGCCACAA	6.3 Mb
274 ~ x	03906-L03066	FANCB gene	Xp22.2	TCTCATCAGAAT-TCTCCTATAAA	114.3 Mb
420 ~ x	00820-L25090	AIFM1 gene	Xq25	TATTGGTCTTGT-GGACAGTAGTTT	25.7 Mb

Length (nt)	SALSA MLPA probe	Gene/Exon <sup>a</sup>	Ligation site	Partial sequence <sup>b</sup> (24 nt adjacent to ligation site)	Distance to next probe
----- Start of PAR2 -----					
355 - +	01156-L00659	VAMP7 gene	Xq28 / Yq12	TGTGGGAAAAGT-GTTTCCATTCTG	98 kb
<i>q</i> -telomere					

<sup>a</sup> See section Exon numbering on page 2 for more information.

<sup>b</sup> Only partial probe sequences are shown. Complete probe sequences are available at [www.mrcholland.com](http://www.mrcholland.com). Please notify us of any mistakes: [info@mrcholland.com](mailto:info@mrcholland.com).

#### Notes:

- The 105 nt chromosome Y-specific control fragment targets the *UTY* gene, located ~11 Mb from the PAR region.
- The 118 nt chromosome Y-specific control fragment targets the *ZFY* gene, located just outside the PAR region at ~470 kb distance from the 456 nt ZBED1 probe. A small signal for the 118 nt fragment may be observed in some female samples.

**CNE** = Conserved Non-coding DNA Element. Locations of the upstream regulatory regions (CNE-2, CNE-3 and CNE-5) are based on Durand et al. (2010) and Benito-Sanz et al. (2012b). Locations of the downstream regulatory regions (CNE2 through CNE9) are based on Benito-Sanz et al. (2012b) and Fukami et al. (2006).

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

- Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

x X-chromosome, outside PAR region. Gives 50% reduced signal in males as compared to females.

∫ A founder *SHOX* point mutation in the Spanish Gypsy population, c.508G>C (p.A170P), results in a strongly reduced probe signal of the 300 nt exon 4 probe (Barca-Tierno et al. 2011).

^ The 481, 318 and 439 nt probes are within the recurrent ~47.5 kb deletion region described by Benito-Sanz et al. (2012b).

+ The VAMP7 probe at 355 nt is located very close to the X and Y *q*-telomere in PAR2.

◇ From product description version 05 onwards the 226 nt probe is considered an exon 6 probe to align with the LRG exon numbering.

# The specificity of this probe relies on a single nucleotide difference between a related gene or pseudogene. As a result, an apparent duplication of only this probe can be the result of a non-significant single nucleotide sequence change in the related gene or pseudogene.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Please note: not all known SNVs are mentioned in the tables above. Single probe aberration(s) must be confirmed by another method.

**Table 3. Examples of final ratios from potential results**

Type of test sample	Probes	RATIO: when compared to		Explanation
		Male ref	Female ref	
Normal FEMALE DNA	All PAR1/PAR2 probes, incl. SHOX AREA, Xp22	1	1	PAR1 & PAR2 regions are located on X and Y: normal females: 2 copies; normal males: 2 copies.
	X-chromosome probes OUTSIDE PAR1/2 region	2	1	Sequences present only on X: normal females: 2 copies; normal males: 1 copy.
	105 & 118 nt Y-chromosome probes outside PAR1 + PAR2	0	-	Sequences present only on Y: normal females: 0 copies; normal males: 1 copy.
Normal MALE DNA	All PAR1/PAR2 probes, incl. SHOX AREA, Xp22	1	1	PAR1 & PAR2 regions are located on X and Y: normal females: 2 copies; normal males: 2 copies.
	X-chromosome probes OUTSIDE PAR1/2 region	1	0.5	Sequences present only on X: normal females: 2 copies; normal males: 1 copy.
	105 & 118 nt Y-chromosome probes outside PAR1 + PAR2	1	∞	Sequences present only on Y: normal females: 0 copies; normal males: 1 copy.
FEMALE DNA with a deletion that includes one or more exons of the SHOX gene	One or more SHOX probes + possibly one or more SHOX AREA and/or Xp22 probes	0.5	0.5	Heterozygous deletion in essential <i>SHOX</i> area in PAR1. PAR1 region is located on X and Y. <b>This typically causes LWD or short stature.</b>
	Deletions extending outside PAR1: one or more of the 6 chromosome X probes	1	0.5	Heterozygous deletion outside PAR1. Sequences present only on X.

Type of test sample	Probes	RATIO: when compared to		Explanation
		Male ref	Female ref	
	(marked x in Table 2) also show a lower copy number			
MALE DNA with a deletion that includes one or more exons of the SHOX gene	One or more SHOX probes + possibly one or more SHOX AREA and/or Xp22 probes	0.5	0.5	Heterozygous deletion in essential SHOX area in PAR1. PAR1 region is located on X and Y. <b>This typically causes LWD or short stature.</b>
	Deletions extending outside PAR1: one or more of the 6 chromosome X probes (marked x in Table 2) also show a lower copy number	0	0	Heterozygous deletion outside PAR1. Sequences present only on X.
FEMALE or MALE DNA with a deletion in the area upstream or downstream of SHOX, but not including SHOX gene probes	One or more probes upstream or downstream of SHOX, indicated with SHOX AREA in table 2	0.5	0.5	Heterozygous deletion. PAR1 region is located on X and Y. <b>Deletions in this region have been associated with LWD &amp; ISS</b> (Benito-Sanz et al. 2005, Benito-Sanz et al. 2012a, Benito-Sanz et al. 2012b, Chen et al. 2009). <b>Not all deletions detected by these probes will result in LWD or ISS!</b>
FEMALE sample with Turner syndrome (45,X)	All PAR1/PAR2 probes, incl. SHOX AREA, Xp22	0.5	0.5	PAR1 & PAR2 regions are located on X and Y. Only one X is present in this case.
	X-chromosome probes OUTSIDE PAR1/2 region	1	0.5	Probe sequences located only on X. Only one X is present in this case.
	105 & 118 nt Y-chromosome probes outside PAR1 + PAR2	0	-	Probe sequences located only on Y. No Y present in this case.
MALE sample with Klinefelter syndrome (47,XXY)	All PAR1/PAR2 probes, incl. SHOX AREA, Xp22.	1.5	1.5	PAR1 & PAR2 regions are located on X and Y. Three instead of the normal two copies present.
	X-chromosome probes OUTSIDE PAR1/2 region	2	1	Probe sequences present only on X. Two Xs are present in this case.
	105 & 118 nt Y-chromosome probes outside PAR1 + PAR2	1	∞	Probe sequences present only on Y. One Y present in this case.
FEMALE sample with Triple X-syndrome (47,XXX).	All PAR1/PAR2 probes, incl. SHOX AREA, Xp22.	1.5	1.5	PAR1 & PAR2 regions are located on X and Y. Three instead of the normal two copies present.
	X-chromosome probes OUTSIDE PAR1/2 region	3	1.5	Probe sequences present only on X. Three Xs are present in this case.
	105 & 118 nt Y-chromosome probes outside PAR1 + PAR2	0	-	Probe sequences located only on Y. No Y present in this case.

## Related SALSA MLPA probemixes

P026 Sotos syndrome	Contains probes for <i>NSD1</i> and <i>NFIX</i> , involved in Sotos syndrome.
P216 Growth Hormone Deficiency mix -1	Contains probes for <i>GH1</i> , <i>PROP1</i> , <i>POU1F1</i> , <i>GHRHR</i> , <i>HESX1</i> , <i>LHX3</i> and <i>LHX4</i> , involved in Growth Hormone Deficiency (GHD).
P217 IGF1R	Contains probes for <i>IGF1R</i> , <i>IGFBP3</i> and <i>IGFALS</i> involved in growth and development.
P262 GHI	Contains probes for <i>IGF1</i> , <i>GHR</i> , <i>JAK2</i> and <i>STAT5B</i> , related to growth hormone insensitivity and short stature.
P329 CRLF2-CSF2RA-IL3RA	Contains probes for the PAR1 genes <i>CRLF2</i> , <i>CSF2RA</i> and <i>IL3RA</i> , linked to B-cell acute lymphoblastic leukaemia (ALL).
P360 Y-Chromosome Microdeletions	Contains probes for the Y chromosomal regions AZFa, AZFb, AZFc, associated with spermatogenetic failure in infertile men.

## References

- Balasubramanian R and Crowley Jr WF (2017). Isolated gonadotropin-releasing hormone (GnRH) deficiency. *GeneReviews [Internet]*. <https://pubmed.ncbi.nlm.nih.gov/20301509/>.
- Barca-Tierno V et al. (2011). Identification of a Gypsy SHOX mutation (p.A170P) in Leri-Weill dyschondrosteosis and Langer mesomelic dysplasia. *Eur J Hum Genet*. 19:1218-25.
- Benito-Sanz S et al. (2005). A novel class of Pseudoautosomal region 1 deletions downstream of SHOX is associated with Leri-Weill dyschondrosteosis. *Am J Hum Genet*. 77:533-44.
- Benito-Sanz S et al. (2006). PAR1 deletions downstream of SHOX are the most frequent defect in a Spanish cohort of Leri-Weill dyschondrosteosis (LWD) probands. *Hum Mutat*. 27:1062.

- Benito-Sanz S et al. (2011a). A non-pathogenic pseudoautosomal region 1 copy number variant downstream of SHOX. *Am J Med Genet A*. 155A:935-7.
- Benito-Sanz S et al. (2011b). Clinical and molecular evaluation of SHOX/PAR1 duplications in Leri-Weill dyschondrosteosis (LWD) and idiopathic short stature (ISS). *J Clin Endocrinol Metab*. 96:E404-12.
- Benito-Sanz S et al. (2012a). Identification of the first PAR1 deletion encompassing upstream SHOX enhancers in a family with idiopathic short stature. *Eur J Hum Genet*. 20:125-7.
- Benito-Sanz S et al. (2012b). Identification of the first recurrent PAR1 deletion in Leri-Weill dyschondrosteosis and idiopathic short stature reveals the presence of a novel SHOX enhancer. *J Med Genet*. 49:442-50.
- Benito-Sanz S et al. (2017). Identification of 15 novel partial SHOX deletions and 13 partial duplications, and a review of the literature reveals intron 3 to be a hotspot region. *J Hum Genet*. 62:229-34.
- Bertorelli R et al. (2007). The homozygous deletion of the 3' enhancer of the SHOX gene causes Langer mesomelic dysplasia. *Clin Genet*. 72:490-1.
- Binder G (2011). Short stature due to SHOX deficiency: genotype, phenotype, and therapy. *Horm Res Paediatr*. 75:81-9.
- Caliebe J et al. (2012). IGF1, IGF1R and SHOX mutation analysis in short children born small for gestational age and short children with normal birth size (idiopathic short stature). *Horm Res Paediatr*. 77:250-60.
- Campos-Barros A et al. (2007). Compound heterozygosity of SHOX-encompassing and downstream PAR1 deletions results in Langer mesomelic dysplasia (LMD). *Am J Med Genet A*. 143A:933-8.
- Chen J et al. (2009). Enhancer deletions of the SHOX gene as a frequent cause of short stature: the essential role of a 250 kb downstream regulatory domain. *J Med Genet*. 46:834-9.
- Dupont C et al. (2007). First reported case of intrachromosomal cryptic inv dup del Xp in a boy with developmental retardation. *Am J Med Genet A*. 143A:1236-43.
- Durand C et al. (2010). Enhancer elements upstream of the SHOX gene are active in the developing limb. *Eur J Hum Genet*. 18:527-32.
- Fukami M et al. (2006). Transactivation function of an approximately 800-bp evolutionarily conserved sequence at the SHOX 3' region: implication for the downstream enhancer. *Am J Hum Genet*. 78:167-70.
- Funari MFA et al. (2019). Evaluation of SHOX defects in the era of next-generation sequencing. *Clin Genet*. 96:261-5.
- Huber C et al. (2006). High incidence of SHOX anomalies in individuals with short stature. *J Med Genet*. 43:735-9.
- Mehta PA and Ebens C (2021). Fanconi Anemia. *GeneReviews [Internet]*. <https://pubmed.ncbi.nlm.nih.gov/20301575/>.
- Mullighan CG et al. (2009). Rearrangement of CRLF2 in B-progenitor- and Down syndrome-associated acute lymphoblastic leukemia. *Nat Genet*. 41:1243-6.
- Nguyen TA et al. (2022). Neuroligins and neurodevelopmental disorders: X-linked genetics. *Front Synaptic Neurosci*. 12:33.
- Russell L et al. (2009). Deregulated expression of cytokine receptor gene, CRLF2, is involved in lymphoid transformation in B-cell precursor acute lymphoblastic leukemia. *Blood*. 114:2688-98.
- Sandbacka M et al. (2011). Evaluation of SHOX copy number variations in patients with Mullerian aplasia. *Orphanet J Rare Dis*. 6:53.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res*. 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat*. 28:205.
- Shears DJ et al. (2002). Pseudodominant inheritance of Langer mesomelic dysplasia caused by a SHOX homeobox missense mutation. *Am J Med Genet*. 110:153-7.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem*. 421:799-801.
- Zinn AR et al. (2002). Complete SHOX deficiency causes Langer mesomelic dysplasia. *Am J Med Genet*. 110:158-63.



## Selected publications using SALSA MLPA Probemix P018 SHOX


- Albuisson J et al. (2012). Clinical utility gene card for: Leri-Weill dyschondrosteosis (LWD) and Langer mesomelic dysplasia (LMD). *Eur J Hum Genet.* 20.
- Auger J et al. (2016). Genotype-Phenotype Relationship in Patients and Relatives with SHOX Region Anomalies in the French Population. *Horm Res Paediatr.* 86:309-18.
- Barroso E et al. (2010). Identification of the first de novo PAR1 deletion downstream of SHOX in an individual diagnosed with Leri-Weill dyschondrosteosis (LWD). *Eur J Med Genet.* 53:204-7.
- Bunyan DF et al. (2013). Diagnostic screening identifies a wide range of mutations involving the SHOX gene, including a common 47.5 kb deletion 160 kb downstream with a variable phenotypic effect. *Am J Med Genet A.* 161A:1329-38.
- Bunyan DJ et al. (2016). Duplications upstream and downstream of SHOX identified as novel causes of Leri-Weill dyschondrosteosis or idiopathic short stature. *Am J Med Genet A.* 170A:949-57.
- Bunyan DJ et al. (2021). Rare dosage abnormalities flanking the SHOX gene. *Egypt J Med Hum Genet.* 22:89.
- Capkova P et al. (2020). Short stature and SHOX (Short stature homeobox) variants-efficacy of screening using various strategies. *PeerJ.* 8:e10236.
- Del Pino M et al. (2019). SHOX deficiency in Argentinean cohort: Long-term auxological follow-up and a family's new mutation. *J Pediatr Genet.* 8:123-32.
- Funari MF et al. (2008). Cryptic intragenic deletion of the SHOX gene in a family with Leri-Weill dyschondrosteosis detected by Multiplex Ligation-Dependent Probe Amplification (MLPA). *Arq Bras Endocrinol Metabol.* 52:1382-7.
- Funari MF et al. (2010). Usefulness of MLPA in the detection of SHOX deletions. *Eur J Med Genet.* 53:234-8.
- Gatta V et al. (2007). Identification and characterization of different SHOX gene deletions in patients with Leri-Weill dyschondrosteosis by MLPA assay. *J Hum Genet.* 52:21-7.
- Gatta V et al. (2014). Spectrum of phenotypic anomalies in four families with deletion of the SHOX enhancer region. *BMC Med Genet.* 15:87.
- Genoni G et al. (2018). Improving clinical diagnosis in SHOX deficiency: the importance of growth velocity. *Pediatr Res.* 83:438-44.
- Hirschfeldova K et al. (2012). SHOX gene defects and selected dysmorphic signs in patients of idiopathic short stature and Leri-Weill dyschondrosteosis. *Gene.* 491:123-7.
- Hirschfeldova K et al. (2017). Detection of SHOX gene aberrations in routine diagnostic practice and evaluation of phenotype scoring form effectiveness. *J Hum Genet.* 62:253-7.
- Kumar A et al. (2020). Pathogenic/likely pathogenic variants in the SHOX, GHR and IGFALS genes among Indian children with idiopathic short stature. *J Pediatr Endocrinol Metab.* 33:79-88.
- Marstrand-Joergensen MR et al. (2017). Prevalence of SHOX haploinsufficiency among short statured children. *Pediatr Res.* 81:335-41.
- Poggi H et al. (2015). A Deletion of More than 800 kb Is the Most Recurrent Mutation in Chilean Patients with SHOX Gene Defects. *Horm Res Paediatr.* 84:254-7.
- Rodriguez FA et al. (2013). Clinical and molecular characterization of Chilean patients with Leri-Weill dyschondrosteosis. *J Pediatr Endocrinol Metab.* 26:729-34.
- Sandoval GT et al. (2014). SHOX gene and conserved noncoding element deletions/duplications in Colombian patients with idiopathic short stature. *Mol Genet Genomic Med.* 2:95-102.
- Stritar J et al. (2021). Detection of del/dup inside SHOX/PAR1 region in children and young adults with idiopathic short stature. *Genes.* 12:1546.
- Stuppia L et al. (2010). Different approaches in the molecular analysis of the SHOX gene dysfunctions. *J Endocrinol Invest.* 33:30-3.






P018 product history	
Version	Modification
G2	One SHOX probe and one SHOX area probe have been adjusted. Minor alteration, no change in sequence detected.
G1	One SHOX area probe has been replaced (CNE4) and four new probes have been included for the upstream SHOX enhancer sequences CNE-2, CNE-3 and CNE-5 and the downstream CNE7 enhancer sequence.
F1	Three new probes near the PAR1 boundary have been included. One probe (GPR143) has been removed. The 88 and 96 nt control fragments have been replaced (QDX2).
E1	Six probes located on chromosome X and six reference probes have been replaced.
D1	One target and six reference probes have been replaced.
C1	Several reference probes replaced. Variable probes in the SHOX downstream region removed. Extra control fragments added.
B1	Many probes outside the SHOX gene have been added.
A1	First release.

Implemented changes in the product description
<p>Version G2-07 – 20 March 2024 (04P)</p> <ul style="list-style-type: none"> <li>-In section Positive control DNA samples, NA04626 was added. The genotype of this sample has been specified in the associated table in that same section.</li> <li>-Morocco has been removed from the list of countries in which the product is IVD-registered.</li> </ul> <p>Version G2-06 – 21 March 2023 (04P)</p> <ul style="list-style-type: none"> <li>- In section Reference samples, clarification added that all reference samples need to be the same sex for correct data analysis.</li> <li>- In section Interpretation of results, expected final ratios for X chromosome probes in males added to the table.</li> <li>- Minor textual and layout changes.</li> </ul> <p>Version G2-05 – 22 February 2022 (04P)</p> <ul style="list-style-type: none"> <li>- Product description adapted to a new template.</li> <li>- Multiple minor textual changes.</li> <li>- Transcript variants section updated according to SHOX LRG_710.</li> <li>- P018 specific notes added about breakpoints in SHOX intron 3, exon 3/intron 3 deletions that are not detected by P018-G2, and note about flanking probes rephrased.</li> <li>- Exon numbering of the SHOX gene has been changed to align with the LRG exon numbering: the 226 nt probe is now considered an exon 6 probe (was intron 6 probe).</li> <li>- Chromosomal band of the 420 nt probe (00820-L25090) updated.</li> <li>- Ligation sites of the probes targeting the SHOX gene updated according to new version of the NM_ reference sequence.</li> <li>- Multiple small changes in Table 2.</li> <li>- Selected publications updated.</li> <li>- UK has been added to the list of countries in Europe that accept the CE mark.</li> </ul> <p>Version G2-04 – 13 November 2020 (02P)</p> <ul style="list-style-type: none"> <li>- Product description rewritten and adapted to a new template.</li> <li>- Various minor textual changes.</li> <li>- Note about a small signal for the ZFY probe in some females added.</li> <li>- Information about slope correction in samples with large deletions added.</li> <li>- Clarification added about the purpose of flanking probes and their relation to other disorders.</li> <li>- Link to mutation database updated.</li> <li>- Note added about copy number variation observed in healthy individuals for a number of probes outside the SHOX region.</li> <li>- SHOX region information in Table 1 updated.</li> <li>- Ligation sites of flanking probes in Table 2 updated.</li> <li>- Small changes of probe lengths in Table 1 and 2 in order to better reflect the true lengths of the amplification products.</li> </ul>

Implemented changes in the product description	
Version G2-03 – 27 April 2020 (04) - Product is now registered for IVD use in Morocco.	
Version G2-02 – 15 April 2019 (04) - Product is now registered for IVD use in Israel.	

More information: <a href="http://www.mrcholland.com">www.mrcholland.com</a> ; <a href="http://www.mrcholland.eu">www.mrcholland.eu</a>	
	MRC Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands
E-mail	<a href="mailto:info@mrcholland.com">info@mrcholland.com</a> (information & technical questions) <a href="mailto:order@mrcholland.com">order@mrcholland.com</a> (orders)
Phone	+31 888 657 200

	EUROPE*  COLOMBIA ISRAEL
	ALL OTHER COUNTRIES

\*comprising EU (candidate) member states and members of the European Free Trade Association (EFTA), and the UK. The product is for RUO in all other European countries.